

Product Warranty

A 6 King Dr Horsham, VIC
P (03) 5381 1404
E sales@oscarfurniture.com.au
W oscarfurniture.com.au

Product Usage

M5 Bariatric Lift-Recline Chair
Applies to: M5 650 & M5 850
Overall weight capacity: 320kg
Footrest weight capacity: 80kg

Proof of Purchase: This Oscar Chair is fitted with a compliance card on the rear bottom rail. All cards have a unique ID number called the "SO" number. Warranty is voided if this card is removed.


Motor/Mechanism: 2 year warranty on all electrical

Foam/Frame: 5-year warranty on material and workmanship, excluding general wear and tear.

Fabric: 2-year warranty on the Versotela fabric ranges when the care labels have been followed
components excluding the battery pack. 1 year warranty on the battery pack.

WARNING:

- The user must ease themselves into the chair in the raised position. dropping into the chair from height may damage the motors and void the warranty.
- Do not place anything over power cords, pass cords through the chair, or drag the chair by the power cord
- Do not allow children to play on this mechanised furniture.



Do not sit on the footrest



Assist footrest by lifting feet



Do not drop into chair when partially raised



Fully raise seat before sitting

Specifications: The power supply voltage is 240 volts AC, 50hz. A transformer reduces the voltage to 24 volts DC. All motors operate at 24V DC

Maintenance: the motors are maintenance free. They are protected by thermal and current overloads and will stop automatically if operated excessively.

Battery Backup: this electrical system is fitted with a 24V Lithium ion battery pack for use in power outages. In the event of a power failure the chair will automatically switch to the emergency backup. The battery pack will automatically recharge when plugged into power.

Cleaning: the hand control has a IPX6 water resistance rating and can be wiped down with a wet cloth but should not be fully submerged in water. Cleaning of fabrics will vary depending on the chosen fabric, for details please check the fabric care card which is sewn into the back of the top head cushion of the chair.

Warranty Terms:

- All warranties apply from the date of sale.
- Where practical the store/distributor will fit warranty replacement parts to minimise repair times for the end user.
- Wherever possible, All warranty claims are to be raised with the store/distributor where the chair was purchased. The store will liaise with the manufacturer to process the claim.

Warranty Exclusions:

- Damage Due to incorrect use of the chair
- General Wear and Tear
- Weight capacity having been exceeded.
- Any modification made to the product by the consumer or a 3rd party without prior written approval from Oscar Furniture to the nature and extent of the modification.
- Minor variations in the finish of the product which does not effect the products value of fitness for its intended purpose.