

Barwon (3-roll back) or Hudson (2-roll back) – B: Spec Sheet

Prepared By	
Client/Organisation name	
Phone	
Address	
Order Date	
Customer PO# (optional)	

Size Modifications:

	Adjustments (mm)				Starting Measurement	Adjustments (mm)						
Seat Depth	-100 <input type="checkbox"/>	-75 <input type="checkbox"/>	-50 <input checked="" type="checkbox"/>	-25 <input type="checkbox"/>	520mm <input type="checkbox"/>	+25 <input type="checkbox"/>	+50 <input type="checkbox"/>					
Back Height		-150 <input type="checkbox"/>	-120 <input type="checkbox"/>	-90 <input type="checkbox"/>	-60 <input type="checkbox"/>	-30 <input type="checkbox"/>	688mm <input type="checkbox"/>	+30 <input type="checkbox"/>	+60 <input type="checkbox"/>	+90 <input type="checkbox"/>	+120 <input type="checkbox"/>	+150 <input type="checkbox"/>
Scapular Height	Linked to Back Height				481mm <input type="checkbox"/>	Linked to Back Height						
Arm Height			-60 <input type="checkbox"/>	-30 <input type="checkbox"/>	190mm <input type="checkbox"/>	+30 <input type="checkbox"/>	+60 <input type="checkbox"/>					
Extended Footrest			-50 <input type="checkbox"/>	-25 <input type="checkbox"/>	N/A <input type="checkbox"/>	+25 <input type="checkbox"/>	+50 <input type="checkbox"/>					
Seat width		-150 <input type="checkbox"/>	-100 <input type="checkbox"/>	-50 <input type="checkbox"/>	540mm <input type="checkbox"/>	+50 <input type="checkbox"/>						
Seat Height		-50 <input type="checkbox"/>	*-30 <input type="checkbox"/>	-20 <input type="checkbox"/>	496mm <input type="checkbox"/>	+25 <input type="checkbox"/>	+45 <input type="checkbox"/>	changed foam, see tips				

Chair Configuration:

*see tips PTO

Mechanism	Single Motor (wall saver) <input type="checkbox"/>	Dual Motor (dual action) <input type="checkbox"/>	Manual <input type="checkbox"/>
	Dual Motor (dual action) must be selected for back height increase of 120 plus		
Castors	Standard Wheel <input type="checkbox"/>	Removed Castors <input type="checkbox"/>	-30 seat height <input type="checkbox"/>
	75mm Wheel <input type="checkbox"/>	B24 Wheel Kit <input type="checkbox"/>	
Battery	Standard 9V Emergency <input type="checkbox"/>	*24V Full Battery Backup <input type="checkbox"/>	(*24V is standard with B24 kit)
Show Wood	Walnut <input type="checkbox"/>	Mahogany <input type="checkbox"/>	Baltic <input type="checkbox"/>
		Teak <input type="checkbox"/>	Remove Show Wood <input type="checkbox"/>
Fabric:			

Customisations:

Memory foam	Seat <input type="checkbox"/>	Leg rest <input type="checkbox"/>	Back rest <input type="checkbox"/>	Arms <input type="checkbox"/>
ROHO Cushion insert	18x16 Low Profile <input type="checkbox"/>	18x16 High Profile <input type="checkbox"/>	16x16 Low Profile <input type="checkbox"/>	16x16 High Profile <input type="checkbox"/>
Controller Side	Standard – (Right) <input type="checkbox"/>			Left <input type="checkbox"/>
Massage Unit	Heated Massage <input type="checkbox"/>			Massage only <input type="checkbox"/>
Glutaral Recess	Standard (No) <input type="checkbox"/>			Yes <input type="checkbox"/>
Leg rest side walls	Standard (No) <input type="checkbox"/>			Yes <input type="checkbox"/>

Accessories:

Headrest (circle size)	A <input type="checkbox"/>	B <input type="checkbox"/>	C <input type="checkbox"/>	D <input type="checkbox"/>	E <input type="checkbox"/>
Head cover	<input type="checkbox"/>				
Arm covers	Fitted <input type="checkbox"/>	Not Fitted <input type="checkbox"/>			
Velcro head & arm covers	<input type="checkbox"/>				

Barwon & Hudson – B: Tips

Seat Height:

Changes in seat height on the B can be achieved through modifications in frame and foam. In the size modification chart under seat depth, changes that require reducing or increasing foam thickness are highlighted in yellow. It is not advised to reduce the seat foam thickness where possible, changing to a Mini size is the preferred option in most cases.

Castors:

- When selecting a B24 wheel kit you cannot change the seat height
- Standard castors are not suitable for moving with the client seated, B24 kit or 75mm castors are required for this. 75mm castor will raise the seat height by 25mm
- B24 Wheel kits include a rear locking brake pedal, push handle and 24V battery backup.
- The no castor option can be done in our factory or on site to make it easier to move into the client’s house. Please indicate the clients preferred option in notes.

General modification notes:

- Gluteal recess: this is a carved-out section in the lower backrest foam for clients with a gluteal region that stops their back from getting adequate contact and support from the chair.
- Battery backup: as standard the single and dual motor chair come with a 9V battery backup capable of lifting the client out of the chair in an emergency. 24V full battery backup will allow the client to use their chair all day without being connected to power.
- Leg rest side walls: a small lip on both sides of the foot rest for clients who have issues with their legs falling off the side of the footrest while reclined.
- Fabric: any fabric can be used on the Hudson chair, contact our office if you require a sample book. Two different fabrics can also be used on one chair i.e. Fusion antibacterial 4-way stretch fabric on contact areas with a hard-wearing Oscar dot on the back and sides.

Additional Notes:

