

## Barwon (3-roll back) or Hudson (2-roll back) – Mini: Spec Sheet

<b>Prepared By</b>	
<b>Client/Organisation name</b>	
<b>Order Date</b>	
<b>Customer PO# (optional)</b>	

### Size Modifications:

	Adjustments (mm)				Starting Measurement	Adjustments (mm)			
<b>Seat Depth</b>	-100 <input type="checkbox"/>	-75 <input type="checkbox"/>	-50 <input type="checkbox"/>	-25 <input type="checkbox"/>	455mm <input type="checkbox"/>	+25 <input type="checkbox"/>	+50 <input type="checkbox"/>	+75 <input type="checkbox"/>	
<b>Back Height</b>		-90 <input type="checkbox"/>	-60 <input type="checkbox"/>	-30 <input type="checkbox"/>	636mm <input type="checkbox"/>	+30 <input type="checkbox"/>	+60 <input type="checkbox"/>	+90 <input type="checkbox"/>	+120 <input type="checkbox"/>
<b>Scapular Height</b>	Linked to Back Height				429mm <input type="checkbox"/>	Linked to Back Height			
<b>Arm Height</b>		-60 <input type="checkbox"/>	-30 <input type="checkbox"/>		199mm <input type="checkbox"/>	+30 <input type="checkbox"/>	+60 <input type="checkbox"/>		
<b>Extended Footrest</b>		-50 <input type="checkbox"/>	-25 <input type="checkbox"/>		N/A <input type="checkbox"/>	+25 <input type="checkbox"/>	+50 <input type="checkbox"/>		
<b>Seat width</b>		-150 <input type="checkbox"/>	-100 <input type="checkbox"/>	-50 <input type="checkbox"/>	540mm <input type="checkbox"/>	+50 <input type="checkbox"/>			
<b>Seat Height</b>		-50 <input type="checkbox"/>	-30 <input type="checkbox"/>	-25 <input type="checkbox"/>	456mm <input type="checkbox"/>	+20 <input type="checkbox"/>	reduced foam, see tips		

### Chair Configuration:

<b>Mechanism</b>	Single Motor (wall saver) <input type="checkbox"/>	Dual Motor (dual action) <input type="checkbox"/>	Manual <input type="checkbox"/>
<b>Castors</b>	Standard Wheel <input type="checkbox"/>	Removed Castors <input type="checkbox"/> -30 seat height	75mm Wheel <input type="checkbox"/>
<b>Battery</b>	Standard 9V Emergency <input type="checkbox"/>	24V Full Battery Backup <input type="checkbox"/>	
<b>Show Wood</b>	Walnut <input type="checkbox"/>	Mahogany <input type="checkbox"/>	Baltic <input type="checkbox"/>
		Teak <input type="checkbox"/>	Remove Show Wood <input type="checkbox"/>
<b>Fabric:</b>			

### Customisations:

<b>Memory foam</b>	Seat <input type="checkbox"/>	Leg rest <input type="checkbox"/>	Back rest <input type="checkbox"/>	Arms <input type="checkbox"/>
<b>ROHO Cushion insert</b>	18x16 (Low Profile) <input type="checkbox"/>		16x16 (Low Profile) <input type="checkbox"/>	
<b>Controller Side</b>	Standard – (Right) <input type="checkbox"/>			Left <input type="checkbox"/>
<b>Massage Unit</b>	Heated Massage <input type="checkbox"/>			Massage only <input type="checkbox"/>
<b>Glutaral Recess</b>	Standard (No) <input type="checkbox"/>			Yes <input type="checkbox"/>
<b>Leg rest side walls</b>	Standard (No) <input type="checkbox"/>			Yes <input type="checkbox"/>

### Accessories:

<b>Headrest (circle size)</b>	A <input type="checkbox"/>	B <input type="checkbox"/>	C <input type="checkbox"/>	D <input type="checkbox"/>	E <input type="checkbox"/>
<b>Head cover</b>	<input type="checkbox"/>				
<b>Arm covers</b>	Fitted <input type="checkbox"/>		Not Fitted <input type="checkbox"/>		
<b>Velcro head &amp; arm covers</b>	<input type="checkbox"/>				

## Barwon & Hudson – Mini: Tips

### Seat Height:

Reduction in seat height on a mini chair is limited due to the size of the mechanism in the chair. Some seat height reduction can be achieved by:

- Use of smaller castors (-25mm)
- Removal of castors completely (-30mm)
- Reduction in seat foam (-20mm) Note: this will void seat foam warranty.
- Max total seat height reduction: -50mm

### Castors:

- The smaller castors will only work effectively on hard surfaces
- Standard and small castors are not suitable for moving with the client seated, 75mm castors are required for this. 75mm castor will raise the seat height by 25mm
- B24 Wheel kit is not available on the Mini.
- The no castor option can be done in our factory or on site to make it easier to move into the client's house. Please indicate the clients preferred option in notes.

### General modification notes:

- Gluteal recess: this is a carved-out section in the lower backrest foam for clients with a gluteal region that stops their back from getting adequate contact and support from the chair.
- Battery backup: as standard the single and dual motor chair come with a 9V battery backup capable of lifting the client out of the chair in an emergency. 24V full battery backup will allow the client to use their chair all day without being connected to power.
- Leg rest side walls: a small lip on both sides of the foot rest for clients who have issues with their legs falling off the side of the footrest while reclined.
- Fabric: any fabric can be used on the Hudson chair, contact our office if you require a sample book. Two different fabrics can also be used on one chair i.e. Fusion antibacterial 4-way stretch fabric on contact areas with a hard-wearing Oscar dot on the back and sides.

### Additional Notes:
